Model:

SE10-45, SE20-45, SE30-45, SE40-45

Brugsvejledning	DK	4
Users manual	UK	9
Mode d' emploi	F	14
Manual de instrucción	ES	19
Manual de Instruções	PT	24

DK Indholdsfortegnelse

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Dixell XR20-30CX (SE-45 Series)

PAR.	DESCRIPTION	RANGE	DEFAULT	U.M.
Set	Set point	LS-US	-45	°C/°F
Hy	diFFerential. Relay compressor tripping differential. The compressor stops on reaching the Setpoint value, and restarts at temperature equal to the Setpoint plus the value of the differential.	0.1...25.0	2,0	°C/°F
LS	Lower SEt. Minimum possible setpoint value.	-55.0...HSE	-45	°C/°F
US	Higher SEt. Maximum possible setpoint value.	LSE...302	-5,0	°C/°F
Ot	CAlibration 1.Calibration 1. Positive or negative temperature value added to the value read by probe 1, based on "ot" parameter settings.	-9.9...9.9	0	°C/°F
P2P	Second probe presence: n=not present, y=present	n/y	N	flag
oE	CAlibration 2.Calibration 2. Positive or negative temperature value added to the value read by probe 2, based on "P2" parameter settings.	-9.9...9.9	0	°C/°F
odS	delay (at) On compressor. Delay time in activating the compressor relay after switch-on of instrument.	0...99	2	min
AC	Anti short cycle delay	0...50	5	Min
CON	On time (compressor). Compressor activation time in the event of faulty probe.	0...99	10	min
COF	OFF time (compressor). Compressor in disabled state time in the event of a faulty probe.	0...99	10	min
CF	Measurement units: °C=Celsius, °F=Fahrenheit	°C/°F	°C	°C/°F
rES	Resolution (only for °C): dE=decimal between -9.9 and 9,9; in=integer	dE-in	dE	flag
Lod	Default display: P1=thermostat probe, P2=evaporator probe, SP=Set point	P1,P2,SP	P1	flag
td	defrost type. Type of defrosting. El = electric defrost. in = reverse cycle defrost (hot gas)	EL/in	In	flag
dtE	Defrost termination temp.	-50...50	6	°C/°F
idF	defrost interval time. Interval between the start of two successive defrosting operations.	0...99	0	hours
MdF	Max length for defrost	0...99	0	min
dFd	Display during defrost: rt=real temp. it=start defrost temp., St=Set point, dF=label Df	Rt/it/SP/dF	St	flag
dAd	Max display delay after defrost	0...255	1	Min
dAF	Defrost delay after fast freezing	0...23	0	hours
ALc	Temperature alarm configuration: rE=relative to SET, Ab=absolute	rE/Ab	rE	Flag
ALU	Higher ALarm. Maximum temperature alarm. Temperature value. which if exceeded in an upward direction triggers the activation of the alarm signal.	AL...99	10	°C/°F
ALL	Lower ALarm. Minimum temperature alarm. Temperature value. which if exceeded in a downward direction, triggers the activation of the alarm signal.	-55...AU	-1	°C/°F
ALd	Temperature alarm delay	0...99	20	min
dAO	Exclusion of temperature alarm at startup.	0...23	1	hour
PbC	Kind of probe	Ptc/ntc	Ntc	flag

Carel Easy Series

	Parameter	Min.	Max.	Def.	UOM
PS	PASSWORD	0	99	22	-
/	PROBE PARAMETERS				
/5	Select °C / °F (0 = °C; 1 = °F)	0	1	0	-
/6	Disable decimal point (1 = disabled)	0	1	1	-
/C1	Probe calibration (OFFSET)	-50.0	50.0	1.8	°C/°F
r	CONTROL PARAMETERS				
St	Setpoint (control temperature)	-50.0	90.0	-45.0	°C/°F
rd	Control differential (hysteresis)	0.0	19.0	3.0	°C/°F
c	COMPRESSOR PARAMETERS				
c0	Comp. and fan start delay after start-up	0	100	0	min
c1	Min. time between successive comp. starts	0	100	3	min
c4	Compressor safety (duty setting)	0	100	15	min
d	DEFROST PARAMETERS				
d0	Type of defrost (0 and 1= defrost by temperature; 2, 3 and 4= defrost by time)	0	4	2	-
dI	Interval between defrosts	0	199	6	h/min
dP	Max. or effective defrost duration	1	199	0	min/s
d4	Defrost when the instrument is switched on (1= activated)	0	1	0	-
d6	Disable temperature display during defrost (1= display disabled)	0	1	1	-
A	ALARM PARAMETERS				
A0	Alarm differential	-20.0	20.0	10.0	°C/°F
AL	Low temperature alarm threshold/deviation	-50.0	250.0	-50.0	°C/°F
AH	High temperature alarm threshold/deviation	-50.0	250.0	250.0	°C/°F
Ad	Low and high temperature alarm delay	0	199	199	min
H	OTHER SETTINGS				
H2	Enable keypad 0= keypad disabled 1= keypad enabled 2= keypad enabled except for ON/OFF function	0	2	1	-
EY	Restore the Default settings	0	1	0	-